
VENO

Gebrauchs- und Pflegeanleitung	2
Medizinische Kompressionsstrümpfe	
Instructions for use and maintenance	4
Medical Pressure Stockings	
Notice d'utilisation et d'entretien	6
Collants médicaux de compression	
Istruzioni per l'uso e manutenzione	8
Calze e collant di contenzione terapeutica	
Instrucciones de uso y mantenimiento	10
Medias de compresión	
Gebruiks- & onderhoudsinstructies	12
Medische steunkousen	
Bruks- og pleieveiledning	14
Medisinske kompresjonsstrømper	

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TIPS FOR PULLING THEM MORE EASILY

1. Pull on your pressure stocking early in the morning after getting up, because then your legs are under "less pressure", i.e. the circumference of the legs is at a minimum. It would be appropriate to wear rubber household gloves when pulling on the stockings, because these gloves give a better grip with the knitted material than bare hands.

2. Turn over the pressure stocking up to the foot part and pull the foot part over the foot. The foot part up to the heel must fit well to the foot. In case of stockings without a toe tip, first of all, it is recommended to put the enclosed foot sliding surface over the foot.

3. Now the pressure stocking can be rolled up the leg carefully until it is pulled on to its full length. The stocking shall not be torn by grasping the edge of the stocking and over-stretching it lengthways.

4. To make the stocking fit well, you should push up the stocking material in "waves" starting from below - doing this, if necessary, in several steps - and finally correcting the fit using the palm of your hand.

WASHING AND MAINTENANCE

Wash your pressure stocking quite often in lukewarm and mild soap water at a temperature of 40°C maximum. After having washed it, rinse the knitted material thoroughly in clear water. Remove the humidity by putting the stocking between two terry towels and pressing it. Your pressure stocking must not be wrung out! For drying, hang the stocking in a shaded area. After drying it, you can iron the knitted material at medium heat with a humid covering cloth.

APPLICATION

The pressure stocking is intended only for medical application in accordance with the medical indication and under medical control. The stocking shall not be worn in the case of the following contra-indications: arterial perfusion disorder, Ulcus cruris, Eczema and wet skin diseases as well as stasis dermatitis. Please consult your doctor if your clinical profile changes to the negative, or if you have unpleasant feelings when using the pressure stockings. Your pressure stocking has been produced in accordance with the quality criteria demanded by the Quality Seal Association. Both compliance with these criteria and the professional delivery of the stocking are guarantees for the pressure stocking's correct pressure in accordance with the prescribed pressure category. Handle your pressure stocking carefully when putting it on or off, and avoid both sharp fingernails and excessive over-stretching. Also ensure that your shoes are in perfect condition.

DURATION OF USE

THUASNE's pressure stockings maintain their medical effect during the entire period of application of six months if they are correctly handled. After that period, your doctor should - if necessary - prescribe a new pressure stocking. If your body size changes significantly due to weight increase or decrease, you should see your consultant doctor.

STORAGE

Store your pressure stocking in such a way that the following damaging influences are avoided: strong sunlight, intense heat, ointments, oils and grease, stain-removing water, petrol, and aggressive washing powder.

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THUASNE DEUTSCHLAND GmbH & Co. KG