

AirCube

CE
0068

INSTRUCTIONS FOR USE

THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Dear customer,

Thank you for choosing Pic Solution's AirCube, an aerosol therapy unit designed and manufactured according to the most up-to-date technologies.

AirCube is a class IIa medical device which converts medicines from liquid or suspension form to aerosol form, so they can be administered directly to the airways.

The nebuliser cup supplied with AirCube allows the medicine to correctly penetrate the respiratory tract, thanks to size of the aerosol particles.

Before using the unit, we highly recommend that you read the short list of warnings contained in this manual, to make sure you have correctly understood how to use it.

ACCESSORIES

A - Adult mask

B - Children's mask

C - Mouthpiece

D - Nose fork

E - Connection tube

F - Nebuliser cup

A, B, C, D, E and F are available as spare parts in kit **REF**
02 00499 000 0000

DEVICE COMPONENTS

G - On/Off switch

H - Compressed air outlet

I - Filter

L - Nebuliser cup support

M - Openings for cooling

N - Power cord

GENERAL WARNINGS

- Before use, carefully read the information in this manual, and keep for future reference.
- Unit designed for aerosol therapy. Follow the instructions of a medical practitioner as to the type of medicine to use, the dose, frequency and duration of inhalations. Any use other than the one for which the unit is intended is considered improper and therefore dangerous; PIKDARE S.r.l. cannot be held liable for any damage caused by incorrect or unreasonable use or if the device is used in electrical systems that do not conform to safety regulations in effect.
- A power outage, a sudden failure or any other adverse condition may cause the unit to stop working, therefore it is recommended that users have an alternative unit or medicine (agreed upon with a physician) that can be used in the event this occurs.
- After removing the device from its packaging, check that it is intact, without any visible damage that may have occurred during transport. If in doubt, do not use the unit, and contact the PIC Service Centre.
- Unroll the power cord completely before connecting it to the mains.
- The incomplete unrolling of the power cord may cause overheating.
- Packaging (bags, box, etc.) may be dangerous and should be kept out of reach of children. Keep the unit out of reach of children, pets and persons with cognitive limitations.
- Before connecting the device, check that your mains voltage corresponds to that shown on the rating label. The rating label is located on the bottom of the unit.
- In case of incompatibility between the unit plug and the electrical socket, use certified adapters if the existing legislation in the individual country allows for it or have an PIKDARE S.r.l. S.p.A. authorised technician or personnel from the PIC service centre replace the plug.
- To prevent overheating and damage to the compressor, switch off the device for at least 40 minutes after every

20 minutes of uninterrupted use.

• WARNING! Using the device before 40 minutes have elapsed may result in motor overheating and cause the safety cut-off devices to activate.

- To ensure correct operation, make sure that the air filter is dry. The filter cover may be removed using a tool with a flat blade. **WARNING! keep the filter cover and the filter out of the reach of children as these are small parts and could represent a swallowing and suffocating hazard.**
- Do not use the device in the presence of nitrous oxide, oxygen or flammable air/anaesthetic mixtures.
- Turn off the appliance and disconnect from the mains after use and before adding more medicine. Do not fill the nebuliser cup beyond the maximum level indicated on it and check that it is assembled with all its components.
- Keep the device and power cord away from hot surfaces.
- Do not use the device while taking a bath or shower, or when in a damp place, or close to bath-tubs, sinks, washbasins, or in any other situation in which there are liquids that may come into contact with the device.
- Never touch the device with wet or damp hands.
- Do not drop or put the device into water or other liquids. Should this occur, unplug it immediately, stop using the appliance and contact PIKDARE S.r.l. S.p.A. or a PIC Service Centre.
- Do not block the openings for cooling air.
- Do not use the device if you feel drowsy.
- If children or persons with physical or mental disabilities use the device, proper supervision by an adult is required. Children must be kept under adult supervision to ensure that they do not play with the device. The device contains small parts that can be swallowed and the power cord poses risk of strangling.
- Do not use connectors or accessories not recommended by the manufacturer.
- Make sure the unit is placed on a level, stable surface when used in order to prevent spillage.
- During use, the surface must be free of any objects that may obstruct the correct flow of air.
- Do not block the holes of the filter cover located on the upper surface of the device with cloth or other items: this could cause the motor to overheat and reduce device performance
- Do not pour liquid on the device, especially on the filter cover: this could lead to its undesired discharge/ nebulisation.
- Before cleaning or performing any maintenance work on the device, disconnect it from the mains by pulling the plug out of the electrical socket.
- If you decide to no longer use the device, disconnect from the mains and cut off the power cord so the device cannot work. Dispose of the device and power cord in compliance with existing laws. We also recommend eliminating any parts that are potential hazards, particularly for children.
- Refer to laws in force for information regarding the disposal of accessories subject to wear. To dispose of the device, refer to Directive 2012/19/EU.
- Do not use the device if, after a fall, it shows signs of damage on any of its parts. If in doubt, do not use the device, and contact the PIC Service Centre.
- In the case of failure or malfunction of the appliance, switch it off by removing the plug from the mains and contact the PIC Service Centre only. The power cord and fuse specified in the "Technical Specifications" must be replaced by PIC Technical Assistance Centre staff only.
- **WARNING! The power cable is the separating element from the electrical mains. When in use, to disconnect the appliance from the mains, switch it off and pull out the plug from the power outlet. Always connect the plug only to a power socket that is easily reachable to facilitate disconnection from the mains. Do not attempt to open, repair or tamper with the appliance.**

PREPARING AND USING THE UNIT

- Take the device and the accessories out of the package.
- Pour the medicine and/or saline solution in the nebuliser tank, according to the doses recommended by your physician, without exceeding the maximum level shown by the mark on the tank; close the nebuliser cup;
- Connect the tube to the nebuliser cup;
- Connect the nebuliser cup to the accessory needed for the therapy: Adult or children's mask, fork, mouthpiece (applied parts) see fig. 3. In case of doubt, ask your physician about which accessory to use.
- Connect the tube to the compressed air outlet and to the nebuliser cup;
- Connect the device to the mains;

- Turn on by pressing the on/off switch and start therapy;
- During the treatment session, the nebuliser cup can be rested on its own support;
- Carry out treatment preferably seated and in a comfortable position;
- When the aerosol jet becomes intermittent, interrupt the treatment for a few seconds to allow the suspended drops of solution to fall on the nebuliser cup walls. Resume the treatment session, which will be finished when no more nebulised solution is released from the nebuliser cup;
- Once the treatment session is completed, switch the device off by unplugging it, remove the accessory used and clean the device and accessories according to the instructions in the Cleaning and Maintenance paragraph.
- The patient and operator are the same person. The steps that the patient/user can do safely are only those specified in the paragraphs "Preparing and using the appliance," "Cleaning and maintenance of the compressor, nebuliser cup and accessories." **Warning!** Do not service and/or repair the appliance while it is being used! Any other maintenance/service activity/operation must be done solely by a PIC Service Centre.

COMPRESSOR CLEANING AND MAINTENANCE

Warning! Do not wet the device or put it into water or other liquids. For cleaning, use a clean, dry cloth only.

CLEANING AND MAINTENANCE OF THE NEBULISER CUP AND ACCESSORIES

Remove and detach all the nebuliser cup components and dip them in boiling water for about 5 minutes.

WARNING! Do not boil the masks.

Rinse the accessories (mask, mouthpiece, fork and tube) under running water. Cold disinfect all parts that come in contact with the patient using denatured alcohol. Carefully dry all accessories and the device before reassembling. Store the device and the accessories in a cool, dry place, away from light and heat. Never use benzene, thinners or other inflammable chemical substances for cleaning.

For increased hygienic safety, never use the same accessories on more than one patient. Purchase a dedicated kit for each user (Aerosol Kit [REF](#)0200499000000).

Clean and check that the condition of the accessories before and after use. Replace them if damaged.

WARNING! Never leave the nose fork unattended as it is a small part that, if swallowed by children, could represent an ingestion and suffocating hazard

CHECKING AND REPLACING FILTERS

Periodically check the condition of the filter. The filter was inserted to protect the compressor. Correct filter maintenance prolongs the life of the aerosol. The frequency with which filter should be replaced depends on the conditions in which the device is used. The filter is located in the upper part of the aerosol.

TO REPLACE THE FILTER:

- Remove the filter cover by lifting it up carefully, using a tool with a flat blade;
- Remove the filter to be replaced from its seat;
- Carefully fit the new filter;
- Reposition the filter cover by pushing it down fully in its seat, making sure that it cannot be removed without the use of a tool.

WARNING! Never leave the removed filter cover and/or filters unattended: These are small parts posing the risk of suffocation if swallowed by children.

TYPE OF DEFECT	CAUSE	SOLUTION
Scarce nebulisation	Clogged nebuliser cup	Remove and detach all the nebuliser cup components and clean them as illustrated in this instruction manual.
Scarce nebulisation	Clogged nebuliser cup	If the issue persists after cleaning, replace the nebuliser cup.
No nebulisation	Nebuliser cup components not assembled correctly	Disassemble and reassemble the nebuliser cup; if the issue persists, replace the nebuliser cup.
No nebulisation	The compressor does not work	Check that the power cord plug is disconnected from the mains and that the on/off switch is on. If the issue persists, contact a PIC Service Centre.
Slow nebulisation	Drug too thick	Dilute the drug as instructed by the manufacturer.

TECHNICAL SPECIFICATIONS

Nominal voltage: 230 V~

Absorbed current: 1A

Maximum flow: 13.5 ± 2 l/min

Max pressure: 2.2 ± 0.6 atm

Conditions for using the device: 20' ON - 40' OFF at 40 °C

MMAD: < 3.7 µm

Nebuliser cup capacity: 10 ml Max ; 2ml Min

Device service life: 1,000 cycles/treatments

Service life of accessories: 1 year

Service life in storage conditions: 1 year

Frequency: 50 Hz

Fuse: T 1 H 250V, 5x20mm

Operational flow: 6 ± 1 l/min

Operational pressure: 0.65 ± 0.2 atm

Weight: 1.2 kg

Breathable fraction: 65.7 % ± 3.4

Noise level at 50 cm*: 68 dB

* Data on noise levels were measured on a new device. Values may vary over time with use.

The information on the unit's capabilities provided by the manufacturer in accordance with EN 13544-1 may not apply to medicines provided in suspension form or those that are highly viscous.

EN 60601-1 Medical electrical equipment – Part 1: general requirements for safety

EN 60601-1-2 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.

- Collateral standard: electromagnetic compatibility.

- Requirements and tests.

EN 13544-1 Respiratory therapy equipment – Part one: Nebulising systems and their components

This section contains information specific to product compliance with the EN 60601-1-2 standard. AirCube is a medical electrical device that requires special precautions regarding electromagnetic compatibility and needs to be installed and commissioned according to the electromagnetic information provided.

Mobile and portable RF communications equipment (mobile phones, transceivers, etc.) may affect the device.

AirCube [REF] 02 038206 000 200 is intended for use in the electromagnetic environment specified below. Customers or users of AirCube [REF] 02 038206 000 200 must ensure that it is used in the recommended environment.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

Emissions test	Complies	Electromagnetic environment - guide
RF emissions CISPR 11	Group 1	AirCube [REF] 02 038206 000 200 uses RF energy only for its internal operation. Therefore, its RF emissions are very low and it is very unlikely that they interfere with nearby electronic devices
RF emissions CISPR 11	Class B	AirCube [REF] 02 038206 000 200 is suitable for use in every room or domestic environment, and in rooms directly connected to a low voltage public network powering buildings used for domestic purposes
Harmonic current emissions IEC 61000-3-2	Class A	
Voltage Fluctuations and Flicker IEC 61000-3-3	Compliant	

Manufacturer's guide and declaration – Electromagnetic Immunity

Immunity test	IEC 60601 test level	Level of compliance	Electromagnetic environment - guide
Electrostatic discharge (ESD) IEC 61000-4-2	In contact ± 6 kV In air ± 8 kV	In contact ± 6 kV In air ± 8 kV	Floors should be wood, concrete or ceramic tile. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power lines	± 2 kV for power lines	The quality of the network voltage must be equivalent to that of a typical commercial or hospital setting.
Surge IEC 61000-4-5	± 1 kV between phases	± 1 kV between phases	
Voltage dips, short interruptions and voltage variations IEC 61000-4-11	<5 % U_1 (>95 % dip in U_1) for 0.5 cycles 40 % U_1 (60 % dip in U_1) for 5 cycles 70 % U_1 (30 % dip in U_1) for 25 cycles <5 % U_1 (>95 % dip in U_1) for 5 s	<5 % U_1 (>95 % dip in U_1) for 0.5 cycles 40 % U_1 (60 % dip in U_1) for 5 cycles 70 % U_1 (30 % dip in U_1) for 25 cycles <5 % U_1 (>95 % dip in U_1) for 5 s	
Power frequency magnetic field (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Electromagnetic fields at network frequency must have levels equivalent to those of a typical commercial or hospital setting.

NB: U_1 is the network AC voltage before applying the test level.

Manufacturer's guide and declaration – Electromagnetic Immunity

<p>Conducted disturbances, induced by radio-frequency fields IEC 61000-4-6</p> <p>Radiated radio-frequency IEC 61000-4-3</p>	<p>3 V eff from 150 kHz to 80 MHz</p> <p>3 V/m from 80 MHz to 2.5 GHz</p>	<p>V_{150} = 3 V eff from 150 kHz to 80 MHz</p> <p>$E_{2.5}$ = 3 V/m from 80 MHz to 2.5 GHz</p>	<p>Portable and mobile RF communication devices must not be used closer to any part of the AirCube [REF] 02 038206 000 200, including the cables, than the separation distance calculated by the equation applicable to the transmitter frequency.</p> <p>Recommended separation distances</p> <p>$d = [3.5 / V_{150}] \sqrt{P}$</p> <p>$d = [3.5 / E_{80}] \sqrt{P}$ from 80 MHz to 800MHz</p> <p>$d = [7 / E_{2.5}] \sqrt{P}$ from 800 MHz to 2.5 GHz</p> <p>Where P is the maximum rated output power of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). The intensities of the fields from fixed RF transmitters, as determined in an electromagnetic study of the site^{a)}, may be less than the level of conformity of each frequency interval^{b)}. Interference may occur near devices marked by the following symbol:</p>
--	---	---	--

NOTE 1 - At 80 MHz and 800 MHz, the distance for the highest frequency range is to be applied.

NOTE 2 - These guidelines may not cover all situations. Electromagnetic radiation is affected by absorption and reflection by structures, objects and people.

^{a)} The field intensities for fixed transmitters like base stations for radio telephones (cellular and cordless) and terrestrial mobile radios, HAM radio devices, AM and FM radio transmitters and TV transmitters can not be theoretically foreseen accurately. In order to assess an electromagnetic environment caused by fixed RF transmitters an electromagnetic study of the site must be considered. If the field intensity measured in the location where the AirCube [REF] 02 038206 000 200 is used exceeds the applicable level of conformity above the normal operation of the AirCube [REF] 02 038206 000 200 must be placed under observation. If abnormal performance is noted additional measures may required such as different orientation or positioning of the AirCube [REF] 02 038206 000 200.

^{b)} The field intensity in 150 kHz to 80 MHz frequency range must be less than $[V_{150}]$ V/m.

More information on distances to be kept between portable and mobile radio communications equipment and AirCube **REF**02 038206 000 200 can be requested from PIKDARE S.r.l. S.p.A., through the contact information provided in this manual. It is nevertheless advisable to keep a safe distance (at least 1 m) between the ultrasound aerosol therapy unit and mobile phones or RF transmitters in order to reduce possible interference.

LEGEND OF SYMBOLS:

Warning! Follow the Instructions for use

I-O On - Off

Class II device

Hz Frequency

IP21

Device protected against penetration of solids and liquids. Device protected against vertically dripping water.

SN Serial Number

Device with Type BF applied part

REF Product identification code

Complies with Directive MDD 93/42/EEC + 2007/47/EC

Manufacturer

Alternate current

Appliance date of production: the second pair of digits in the lot number identifies the year, the third pair identifies the month. (LOT 1016100001; 16 = 2016; 10 = October)

CONDITIONS FOR STORAGE AND TRANSPORTATION:

Relative humidity: 15% / 85% RH

Temperature: -25°C / 70°C

Store in a cool, dry place.

ENVIRONMENTAL CONDITIONS DURING USE:

Relative humidity: 15% / 85% RH

Temperature: 5°C / 40°C
Atmospheric pressure from 700hPa to 1060hPa

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

WARRANTY

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions for use.

The warranty shall not therefore apply in the case of damages caused by improper use, wear or accidental events. For the duration of warranty on conformity defects please refer to the specific provisions of national laws applicable in the country of purchase, where provided.

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