# VENOFLEX®

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Conserver cette notice.

Dear patient,


Your doctor has prescribed you VENOFLEX compression socks, stockings or tights. These are specially designed to relieve the venous problems that you are suffering from, both preventing and treating them.

VENOFLEX compression socks, stockings or tights are at the cutting edge of technology and allow you to combine effectiveness, comfort and aesthetic properties. Regularly wearing them will guarantee successful treatment.

### INSTRUCTIONS FOR PUTTING THEM ON

The socks, stockings or tights should preferably be put on in the morning, whenever you get up, since your legs are not swollen then. Avoid direct contact of your socks, stockings or tights with creams or ointments.

#### Compression classes 1 and 2

*It is recommended that you wear rubber gloves to make it easier to slip on your socks, stockings or tights. Be careful not to damage them with your nails.*

*Roll your socks, stockings or tights back on themselves, right to the foot. Slip on, then check that the heel is in the right position.*


*Using both hands, stretch your socks, stockings or tights widthwise to make it easier to pull them over the ankle and carefully unroll them up your leg. Avoid pulling on the fabric any more than necessary.*

*Adjust your stockings or tights properly at the knee to make sure that there are no wrinkles along the leg.*

#### Compression classes 3 and 4

*It is recommended that you wear rubber gloves to make it easier to slip on your socks, stockings or tights. Be careful not to damage them with your nails.*

*Turn the socks, stockings or tights back on themselves, right to the heel.*

*Pull the turned-back socks, stockings or tights over your foot up to the heel then pull over the heel and onto your leg. Check that the heel is in the right position.*

*Adjust your stockings or tights properly at the knee to make sure that there are no wrinkles along the leg.*

### WASHING INSTRUCTIONS

1. Wash daily using a detergent for delicate fabrics. Do not use fabric-softener.
2. Rinse thoroughly with warm water.
3. Squeeze out excess water without twisting and do not use a tumble-drier or direct heat sources (radiator, direct sunlight, etc.) to dry.
4. Do not iron.

### INDICATIONS

- Evening oedema
- Varicosity
- Telangiectasia
- Simple varicose veins
- Complicated varicose veins or with a previous history
- Following sclerotherapy
- Following surgery
- Leg ulcers
- Skin conditions
- Venous insufficiency during pregnancy
- Prevention of deep vein thrombosis
- Deep vein thrombosis (established accident)
- Post-thrombotic disease
- Lymphedema

### DURATION OF USE

If you follow the instructions for putting on and washing your socks, stockings or tights, they will retain their therapeutic efficacy for a long time. However, given that they are worn regularly and washed frequently, it is recommended that you renew them every 6 months, after consulting your doctor.

### VENOFLEX MICROFIBRES

These are new-generation compression stockings or tights, made with MICROFIBRES. In most cases, it is possible to repair a pulled thread by exerting traction around the stitch.

### STOCKINGS WITH ANTI-SLIP (HOLD-UP STOCKINGS)

To retain the anti-slip properties, make sure that you correctly follow the washing instructions above. Positioning the stockings at the very top of the leg helps to ensure good hold.

When putting them on, do not pull the anti-slip part to avoid any damage.

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